NSK

OPERATION MANUAL

Please read this OPERATION MANUAL carefully before use, and file for future reference.

Contra Angle Handpiece


FX15 / FX15m
FX22 / FX23 / FX23m
FX25 / FX25m
FX75 / FX75m

FX15 / FX15m / FX22 / FX23 / FX23m
FX25 / FX25m / FX75 / FX75m


CE 0197

MADE IN JAPAN


OM-C0438E 002


1


2


3


4


5


6


7


8


11


12


13


14

15

16

17

1. User and Indications for Use

User: Qualified Professionals

Indications for Use:

FX15 / FX15m / FX22 / FX23 / FX23m / FX25 / FX25m

This medical device is intended for transmitting rotation of the power source at the direct drive ratio or at a different gear ratio, thereby running the instrument such as a bur or a reamer to cut and polish natural or artificial teeth during dental treatment.

FX75 / FX75m

This medical device is intended for transmitting rotation of the power source at a reduction gear ratio, thereby running the rotary endodontic instrument, etc. to perform endodontic treatments such as root canal enlargement.

2. Precautions for handling and operation

- Please read these precautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

Class	Degree of Risk
 WARNING	Hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.
 CAUTION	Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.
NOTICE	General product specification information highlighted to avoid product malfunction and performance reduction.

 **WARNING**

- Immediately after a treatment, the handpiece should be cleaned, lubricated and sterilized. Failure to properly maintain the handpiece may lead to overheating, causing burn injuries or product failure. (Refer to "8. Post-use Maintenance")

- If blood infiltrates inside a handpiece, an automatic handpiece cleaning and lubrication system may not totally clean the internal handpiece components, which may lead to internal blood coagulation. Coagulated blood inside the handpiece may cause handpiece failure and overheating, causing burn injuries. To avoid the risk, NSK recommends to immediately lubricate those handpieces with PANA SPRAY Plus / PANA SPRAY.
- Depressing the Push Button while the handpiece is in rotation may lead to overheating, causing burn injuries or product failure. Avoid the push button to contact with any oral tissue. (Except for FX22)

- Supply coolant water and coolant air while using the handpiece. No supplying the coolant water and coolant air may lead to overheating, causing burn injuries or product failure.
- Do not use the handpiece for polishing purposes. Polishing paste could enter the handpiece, causing the failure of the Push Button or burn injuries due to the overheating of the handpiece head.
- Keep any debris or other foreign materials away from inside the gear or the handpiece. Foreign materials remaining inside may lead to overheating, causing burn injuries or other accidents.

- Ball bearings are wear items and may require replacement. Be sure to inspect the handpiece before use (Refer to "7. Checking the Handpiece Before Each Use"). If abnormal vibration or noise are found, the bearings or other internal parts may be worn out or damaged. Damaged parts may lead to handpiece overheating causing burn injuries. If overheating is suspected, stop using the handpiece immediately and contact your Authorized NSK Dealer.

CAUTION

- When operating the product always consider the safety of the patient.
- Users are responsible for the operational control, maintenance and continual inspection of this product.
- Prior to clinical use, inspect the handpiece. Check for vibration, noise and overheating. If any abnormalities are found, stop using the handpiece immediately and contact your Authorized NSK Dealer. (Refer to "7. Checking the Handpiece Before Each Use")
- Do not disassemble or alter the handpiece except as recommended by NSK in this Operation Manual.
- Do not allow any impact on to the product. Do not drop the product.
- Operators and all others in the area must wear eye protection and a mask when operating this handpiece.
- Do not use burs/files with problems listed below as the bur/file may break, seize up or disengage from the chuck.

- Bent, deformed, worn, rusted, broken, deficient bur/file.
- Bur/File which is cracked on the edge or axis.
- Non-ISO(EN ISO) standard, or tampered bur/file.
- Do not use burs/files other than specified in "11. Specifications" as such burs/files may be accidentally released or break during rotation.
- Always keep the bur/file shank clean. Dirt or debris in the chuck could cause poor bur/file concentricity or low chuck retention force.
- Always insert the bur/file all the way into the chuck. If insertion is insufficient, premature failure of the bearings or accidental release of the bur/file may occur.
- Do not exceed Max. Bur Length shown in "11. Specifications".
- Always follow the instructions provided by the bur/file manufacturer.
- Do not exceed the bur/file speed recommended by the bur/file manufacturer.

- Should the handpiece function abnormally during use, stop using the handpiece immediately and contact your Authorized NSK Dealer.
- Do not use the following fluids to wipe, immerse or clean the product; strong/super acid water, strong acid/alkaline chemicals, chlorine-containing solutions, solvents such as benzine or thinner.
- The handpiece is delivered in a non-sterile condition and must be autoclaved prior to use.
- Perform regular function and maintenance checks. (Refer to "10. Periodical Maintenance Checks")
- If the handpiece has not been used for a long period, rotate the handpiece and check for noise, vibration and overheating before use.
- To avoid clinical downtime it is recommended that a spare be kept on hand in case of a breakdown during treatment.

- The operation of the handpiece is permitted only on dental units which correspond to the standards IEC 60601-1 (EN 60601-1) and IEC 60601-1-2 (EN 60601-1-2).
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

3. Accessory List

No.	Part Name	Quantity	Remarks
1	FX Wrench	1	Included in the package for FX22 / FX23 / FX23m / FX25 / FX25m
2	Cap Wrench	1	Included in the package for FX22
3	E-Type Spray Nozzle	1	-

4. Connecting and Disconnecting the Handpiece from the Motor

4-1 Connecting

- 1) Insert the handpiece direct into the motor (Fig. 1).
- 2) Confirm that the handpiece is securely connected to the motor.

4-2 Disconnecting

Hold the motor and the handpiece, then pull apart.

 **CAUTION**

- Do not connect or disconnect the handpiece until the motor has completely stopped.
- Connect ONLY to E type motors (ISO 3964 (EN ISO 3964)).
- Do not exceed the Max. Rotation Speed (Motor) shown in "11. Specifications".

5. Inserting and Removing the Bur/File

5-1 Inserting the Bur/File

FX15 / FX15m / FX23 / FX23m / FX25 / FX25m / FX75 / FX75m (Fig. 2)

- 1) Insert the bur/file until it is correctly secured in place.
- 2) Depress the Push Button and insert the bur/file into the chuck until the bur/file "notch" mechanism engages. Release the button.
- 3) Confirm that the bur/file is secure by gently pulling and pushing the bur/file without depressing the Push Button.

 **CAUTION**

- Only FX75 / FX75m can use file.

FX22 (Fig. 3, 4)

- 1) Turn the Latch Plate to open the chuck.
- 2) Insert the bur into the chuck until the bur securely engages. Swing the Latch Plate back to closed position.
- 3) Ensure that the bur is secure by gently pulling and pushing the bur.

5-2 Removing the Bur/File**FX15 / FX15m / FX23 / FX23m / FX25 / FX25m / FX75 / FX75m** (Fig. 5)

Depress the Push Button to open the chuck (❶) and remove the bur/file (❷).

FX22 (Fig. 3)

Turn the Latch Plate to open position and pull the bur out.

 **CAUTION**

- Do not use burs/files with problems listed below as the bur/file may break, seize up or disengage from the chuck.
 - Bent, deformed, worn, rusted, broken, deficient bur/file.
 - Bur/File which is cracked on the edge or axis.
 - Non-ISO(EN ISO) standard, or tampered bur/file.
- Do not use burs/files other than specified in "11. Specifications" as such burs/files may be accidentally released or break during rotation.

- Always keep the bur/file shank clean. Dirt or debris in the chuck could cause poor bur/file concentricity or low chuck retention force.
- Always insert the bur/file all the way into the chuck. If insertion is insufficient, premature failure of the bearings or accidental release of the bur/file may occur.
- Do not exceed Max. Bur Length shown in "11. Specifications".
- Always follow the instructions provided by the bur/file manufacturer.
- Do not exceed the bur/file speed recommended by the bur/file manufacturer.
- Do not insert or remove the bur/file until the rotation has completely stopped.
- Do not apply excess pressure to the bur/file as it may break or bend or become difficult to remove.

6. Connection of the Spray Nozzle (FX22 / FX23 / FX23m / FX25 / FX25m)

- 1) Attach the Spray Nozzle (Optional) to the handpiece head (Fig. 6).
 - 2) Cut the Irrigation Tube to the desired length then connect to the Spray Tube Nozzle on the motor (Fig. 6).
- *Refer to "14. Option Parts List" to identify the correct cartridge.

7. Checking the Handpiece Before Each Use

Follow the check procedure below before use. If any abnormalities are found, stop using the handpiece immediately and contact your Authorized NSK Dealer.

- 1) Check the Head Cap or Cap (FX22) is firmly tightened.
- 2) Insert the bur/file. (Refer to "5. Inserting and Removing the Bur/File")
- 3) Rotate the handpiece for about one minute at the Max. rotation speed of the attached motor. During rotation, check for abnormalities such as abnormal rotation, vibration, noise.
- 4) After the handpiece rotation has completely stopped, touch the handpiece head to confirm the head is NOT heating abnormally (Fig. 7).

CAUTION

- To avoid injury, keep your hands away from the bur during rotation.

8. Post-use Maintenance

After each patient maintain the product as follows.

WARNING

- Immediately after a treatment, the handpiece should be cleaned, lubricated and sterilized. Failure to properly maintain the handpiece may lead to overheating, causing burn injuries or product failure.

CAUTION

- Do not use the following fluids to wipe, immerse or clean the product; strong/super acid water, strong acid/alkaline chemicals, chlorine-containing solutions, solvents such as benzine or thinner.

8-1 Preparation

- 1) Wear a eye protection, a mask and gloves to avoid infection.
- 2) Remove the bur/file.
- 3) Disconnect the handpiece from the motor.

8-2 Cleaning

■ Manual Cleaning (External Surface)

- 1) Remove dirt and debris on the handpiece surface using a brush under running water (at 38°C or less). Do not use a wire brush.
- 2) Wipe off the moisture on the surface using absorbent cloth.

CAUTION

- If water comes inside the handpiece, blow air into the rear of the handpiece using an air syringe and remove the remaining water.

■ Automatic Cleaning (External and Internal Surfaces) (Only for FX15m / FX23m / FX25m / FX75m)

Use a thermo-disinfector in accordance with ISO 15883 (EN ISO 15883).

Refer to the thermo-disinfector manual.

CAUTION

- After washing with thermo-disinfector and prior to lubrication, dry the product until all internal moisture is thoroughly removed. Thermo-disinfector moisture remaining inside the product could reduce the effect of lubrication and could cause corrosion inside of the product.

8-3 Disinfection

■ Manual Disinfection (External Surface)

Wipe the surface with cloth soaked in disinfectant solution such as alcohol.

■ Automatic Disinfection (External and Internal Surfaces) (Only for FX15m / FX23m / FX25m / FX75m)

Use a thermo-disinfector in accordance with ISO 15883 (EN ISO 15883).

Refer to the thermo-disinfector manual.

8-4 Lubrication

■ NSK PANA SPRAY Plus / PANA SPRAY

Apply NSK PANA SPRAY Plus / PANA SPRAY every time after each use and/or before autoclaving.

- 1) Attach the Spray Nozzle into the Spray Port nozzle on the can.
- 2) Insert the Spray Nozzle into the rear of the handpiece. Hold the handpiece and spray for approximately 2-3 seconds. Apply lubricant until it expels from the handpiece head for at least 2 seconds (Fig. 8).

⚠ CAUTION

- Hold the spray can upright.
- Be sure to hold the handpiece firmly to prevent slipping out of hand by the spray pressure when lubricating.
- Supply lubricant until it comes out of the handpiece head.
- If you wish to expel excessive oil from inside of the handpiece, rotate the handpiece for approx. 15 seconds without a bur/file. During the rotation, do not depress the push button inadvertently, especially when using absorbent cloth to prevent oil from scattering. Depressing the push button during rotation may cause a failure to attach a bur/file. (Except for FX22)


Fig. 9


Fig. 10

NOTICE

- NSK recommends the use of "Spray Mist Absorber" (Y900084) to prevent oil mist expelling out of the handpiece head.

■ NSK automatic handpiece cleaning and lubrication system

When using NSK automatic handpiece cleaning and lubrication system, refer to the system instructions.

WARNING

- If blood infiltrates inside a handpiece, an automatic handpiece cleaning and lubrication system may not totally clean the internal handpiece components, which may lead to internal blood coagulation. Coagulated blood inside the handpiece may cause handpiece failure and overheating, causing burn injuries. To avoid the risk, NSK recommends to immediately lubricate those handpieces with PANA SPRAY Plus / PANA SPRAY.

8-5 Sterilization

Sterilize the product by autoclave sterilization. After each patient sterilize as noted below.

1) Insert into an autoclave pouch. Seal the pouch.

2) Autoclavable under the conditions below.

Autoclave for more than 20min. at 121°C, or 15min. at 132°C, or 3min. at 134°C.

3) The handpiece should remain in the autoclave pouch until required for use.

CAUTION

- Be sure to clean and lubricate the handpiece prior to autoclaving. If blood remains on the internal surface it can become clotted and cause product failure.
- Do not autoclave the product with other instruments even when it is in a pouch. This is to prevent possible discoloration and damage to the product from chemical residue on other instruments.
- Do not heat or cool the product too quickly. Rapid change in temperature could cause damage to the product.
- To avoid product failure, do not use a sterilizer whose temperature exceeds 138°C during use including the dry cycle. The chamber temperature of some type of sterilizers may exceed 138°C. Contact the sterilizer manufacturer for detailed information.

- Keep the product in suitable atmospheric pressure, temperature, humidity, ventilation, and sunlight. The air should be free from dust, salt and sulphur.
- Do not touch the product immediately after autoclaving as it will be very hot and must remain in a sterile condition.
- Autoclave sterilization is recommended for the product. The validity of other sterilization methods is not confirmed.

NOTICE

- NSK recommends Class B sterilizers as stated in EN 13060.

9. Replacing the Cartridge (FX22 / FX23 / FX23m / FX25 / FX25m)

- 1) Detach the shank by turning it counter clockwise and remove the head from the shank (Fig.11).

CAUTION

- After removing, do not hold the shank upside down as it could cause the internal parts to drop.

- 2) Insert the wrench firmly into the rear of the handpiece head until it stops, as shown in Fig. A-1.
- 3) Tightly hold the FX wrench as illustrated in Fig. B then pull the drive shaft straight out in the direction of the arrow (Fig. B).


Fig. A-1


Fig. A-2


Fig. B

CAUTION

- Be sure to insert the FX wrench straight into the handpiece head. Do not insert the FX wrench at an angle as shown in Fig. A-2.
- When removing the drive shaft, tightly hold the FX wrench as directed in Fig. B.

FX23 / FX23m / FX25 / FX25m

- 4) Locate the correct wrench tool for the head cap then turn the wrench counter clockwise to loosen the cap. Remove the cap (Fig.12).
 - 5) Use the bur to gently push the entire cartridge out from the head (Fig. 13).
 - 6) Clean the head interior with NSK PANA SPRAY Plus / PANA SPRAY.
 - 7) Wipe the NSK PANA SPRAY Plus / PANA SPRAY oil over the head interior.
 - 8) Insert the new cartridge into the head by aligning the pin on the cartridge with the slot on the head (Fig. 14).
 - 9) Tighten the head cap with the correct wrench.
 - 10) Insert the drive shaft firmly into the handpiece head until it clicks. If the click is not heard, insert it again.
 - 11) Reattach the shank.
- *Refer to "15. Spare Parts List" to identify the correct cartridge.

FX22

- 4) Locate the correct Cap Wrench tool for the cap then turn the cap wrench counter clockwise to loosen the cap. Remove the cap (Fig. 15).
 - 5) Use the bur to gently push the entire cartridge out from the head (Fig. 16).
 - 6) Clean the head interior with NSK PANA SPRAY Plus / PANA SPRAY.
 - 7) Wipe the NSK PANA SPRAY Plus / PANA SPRAY oil over the head interior.
 - 8) Insert the new cartridge into the head by aligning the pin on the cartridge with the slot on the head (Fig. 17).
 - 9) Tighten the cap with the correct cap wrench.
 - 10) Insert the drive shaft firmly into the handpiece head until it clicks. If the click is not heard, insert it again.
 - 11) Reattach the shank.
- *Refer to "15. Spare Parts List" to identify the correct cartridge.

CAUTION

- Before turning the head cap, make sure the wrench is correctly aligned with the cap. Do not turn the wrench at an angle as abrasion of the wrench may occur.
- Secure the shank firmly to the handpiece head so that there is no gap between them. A loose connection could make the head unstable, resulting in an accident during use.
- Use only a genuine NSK cartridge. If a non genuine cartridge is used NSK cannot guarantee performance and the handpiece warranty would become invalid.
- Always first finger tighten the head cap then secure firmly with the wrench.
- NSK never recommends the disassembly and repair of any NSK cartridge. There is NO EXCEPTION. In such a case a handpiece may perform abnormally (abnormal noise or abnormal vibration). Damage, failure or accidents are outside of our guarantee.

10. Periodical Maintenance Checks

Every 3 months perform periodical maintenance checks, referring to the check sheet below. If any abnormalities are found, contact your Authorized NSK Dealer.

Points to check	Details
Head Cap or Cap (FX22) is loose	Tighten firmly using the correct head cap wrench or cap wrench.
Rotation	Rotate the handpiece and check for abnormalities such as abnormal rotation, vibration, noise, and overheating.

11. Specifications

Model	FX15	FX15m	FX22	FX23	FX23m	FX25	FX25m	FX75	FX75m
Max. Rotation Speed (Motor)	40,000 min ⁻¹		30,000 min ⁻¹			40,000 min ⁻¹		20,000 min ⁻¹	
Max. Rotation Speed (Handpiece)	10,000 min ⁻¹		30,000 min ⁻¹			40,000 min ⁻¹		1,250 min ⁻¹	
Gear Ratio	4:1 Reduction		1:1 Direct Drive				16:1 Reduction		
Bur Type	ISO 1797-1 (EN ISO 1797-1) Type1 Ø2.35mm CA Bur							ISO 1797-1 (EN ISO 1797-1) Type1 Ø2.35mm Engine File	
Chuck Type	Push Button Type		Latch Type	Push Button Type					
Chuck Length	12.7mm		13.2mm	12.7mm			9.4mm		
Max. Bur Length	22.5mm							-	
Max. Working Part Diameter	Ø4.0mm							-	
File Motion	-							360°Rotation	
Thermo Disinfector	-	X	-	X	-	X	-	X	
Use Environment	Temperature: 10 - 40°C, Humidity: 30 - 75% (No Condensation)								
Transportation and Store Environment	Temperature -10 - 50°C, Humidity: 10 - 85%, Pressure: 500 - 1,060hPa								

* Available models differ according to countries/regions.

12. Symbol

This product is Autoclavable up to Max.135°C.

This product can be washed via thermo disinfectior.

0197

Conforms to CE European Directive of "Medical device directive 93/42/EEC."

Manufacturer.

Authorized representative in the European community.

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

13. Warranty

NSK products are warranted against manufacturing errors and defects in materials. NSK reserves the right to analyze and determine the cause of any problem. Warranty is voided should the product be not used correctly or for the intended purpose or has been tampered with by unqualified personnel or has had non NSK parts installed. Replacement parts are available for seven years beyond discontinuation of the model.

14. Option Parts List

Model	Order Code	Compatible Product
FX Spray Nozzle	C1052702	FX15 / FX15m / FX23 / FX23m / FX25 / FX25m
Spray Nozzle	C032701	FX22

15. Spare Parts List

Model	Order Code	Compatible Product	Remarks
FX23-03	C1056700	FX23 / FX23m	Cartridge
FX25-03	C1052700	FX25 / FX25m	Cartridge
NAC-03	C030011	FX22	Cartridge
FX25-04	C1052500	FX23 / FX23m / FX25 / FX25m	Head Cap
FX25-02	10000946	FX23 / FX23m / FX25 / FX25m	Drive Shaft
FX22-02	10000996	FX22	Drive Shaft
FX Wrench	C1052701	FX23 / FX23m / FX25 / FX25m	-
Cap Wrench	C071751	FX22	-
E-Type Spray Nozzle	Z019090	-	-

16. Disposing product

In order to avoid the health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose the product for you.

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NSK Europe GmbH
NSK France SAS
NSK United Kingdom Ltd
NSK America Corp
NSK Oceania Pty Ltd